

WATER DISTILLATION / DZ SERIES

Stainless-Steel Double Water Distiller

DZ5C / DZ10C / DZ20C



FEATURES

- Stainless-steel construction with two distillation vessels.
- First and second distillation stages with two heating control systems.
- Seamless stainless-steel condenser tubing supports effective cooling.
- Heating stops automatically when the water level is low and resumes when the normal level is restored.
- Copper immersion heating tubes provide efficient heat transfer.

TECHNICAL PARAMETERS

Model	DZ5C	DZ10C	DZ20C
Power supply	380 V / 50 Hz	380 V / 50 Hz	380 V / 50 Hz
Heating power	10 kW	16 kW	32 kW
Water output	≥5 L/h	≥10 L/h	≥20 L/h
Wooden packing size (cm)	45×68×84	46×77×95	56×86×120
Gross weight	28 kg	38 kg	59 kg
Net weight	13 kg	18 kg	30 kg

Product images and specifications are taken from the supplied BIOBASE catalogue, page 18.